

# GETINGE QUADRO BENEFITS

## Sterilizes more in less space and time

- Standard B-process in just 25 minutes, incl. drying!
- More instruments back in circulation faster
- Need of less expensive instruments - *cuts your costs!*



## Process release and traceability

- Process and load release directly on the display
- Cost effective traceability without any special programs
- Label printer connected to sterilizer
- Possible process storage via Ethernet
- Use the clinic computer for process storage  
- *no extra computer needed!*



## Superior loading capacity

- Quadric chamber for maximum loading capacity, e.g. five wrapped dental cassettes
- Fits trays, containers etc. with max depth 313 mm
- Better loading capacity - *increased productivity!*



## Small footprint

- Small footprint and weight of ~50 kg
- Easy to install and fits in limited areas
- *Very high capacity to footprint!*



ECLIPSE DENTAL ENGINEERING  
Westview Court  
10 -14 Kent Road,  
Dartford, Kent DA1 2AJ

Tel: 01322 293333  
Fax: 01322 421158  
Email: [info@eclipse-dental.com](mailto:info@eclipse-dental.com)  
Website: [www.eclipse-dental.com](http://www.eclipse-dental.com)

