

Silent PRIME versions without air dryer

MINI BOX 24/7 PRIME S











Code: PR-MB100-024-15





Description: Oil free compressor, super silenced, FAD at 5 bar 85 l/min, tank capacity 24 l., without air dryer.

Complete with MINI BOX acoustic cabinet that can be installed also separately and in working environments.

MINI BOX 2016 acoustic cabinet provides an optimal noise reduction (50 dB/1M)

User friendly, compact and light, it can be placed everywhere.

Model									
	Dental chairs	Air dryer	Air Capacity m3/h l/min cfm	FAD @ 5bar m3/h l/min cfm	Motor HP KW	Volt. V	Freq. Hz	Receiver l GAL	
MINI BOX 24/7 PRIME S	1	no	- 120 4,2	- 85 3,0	1 0,75	230	50	24	6,34

Model		Duty cycle					
	Noise Level dB(A)		bar	PSI	Size mm	Net weight Kg lb	
MINI BOX 24/7 PRIME S	50	70%	7	101	530x500x800	65	143

