

Long sleeved apron

Waterproof and thumb looped for maximum protection



These long sleeved thumb looped aprons are ideal for preventing and isolating dust, particles, alcohol, blood, bacteria and viruses.

Being waterproof, they are widely used in hospitals, medicine checks, laboratories, hotels, housework, painting, coating, printing, dying, food processing etc

Comfortable and generous in size, the apron easily ties at the back and encompasses thumb loops on the long sleeves to maximise protection of the arms, as well as the body.

- ✔ Waterproof
- ✔ Long sleeved
- ✔ Thumb loops
- ✔ Open back gown
- ✔ Material: CPE
- ✔ Colour: Blue
- ✔ Size: 113 x 107 x 193cm



The aprons are delivered in 5 pieces per poly bag with 300 pieces in a carton. Carton size: 41 x 50 x 43cm with a gross weight of 16.5kg.

RoHS Directive (EU) 2015/863 amending Annex II to Directive 2011/65/EU Test Method: With reference to IEC 62321-4:2013+A1:2017, IEC62321-5:2013, IEC62321-7-2:2017, IEC62321-6:2015 and IEC62321-8:2017, analysed by ICP-OES, UV-Vis and GC-MS