

dental innovations

Filled with cutting edge products and promotions for dentists

IN THIS ISSUE:

Ancar dental treatment centres
Space saving air compressor
and suction pump
Brighten up your surgery

CrossTex DentuaPure® DP365B
One bur does it all
EclipseCare
Celebrating 12 years

CASE STUDY



A State-of-the-art transformation

Dr Jonathan Hall has built a reputation for excellent standards in dental care and the sheer number of patients who have travelled internationally to experience his skill and expertise testify to this fact. He has been appointed as the personal dentist to a number of leading members of state across the world and the preferred dentist for a number of foreign missions based in London and receives patient referrals worldwide.

Dr Hall discussed his 'One Harley Street' practice's recent refurbishment with The Dentist magazine (November Edition).

"We had a couple of 20-year old Sirona units at the practice, and although they were functional and reliable, they were not modern in appearance. So we decided to update one of the units and refurbish the surgery to make the workflow easier, and have a more up to date look."

We were recommended by Bjorn Friis

of Heka Dental as a reputable dealer who could supply, install and service the treatment centre. After interviewing other dentists whom we had done work for, he was satisfied that we would be able to deliver the type of solution he wanted.

We created a detailed plan based on Jonathan's initial ideas and drawings and worked to a tight schedule. The raised floor was removed to eliminate a step in to the surgery and services were moved below the original floor. The surgery was completely refurbished with the air-conditioning unit relocated.

We supplied and installed all dental cabinetry and equipment including a Heka Dental UNIC Treatment Centre.

The surgery has been transformed into a modern and luxurious space be-fitting for the clientele experience and the professionals who work here.



“ We are very pleased with the transformation of the space, and all the planning and skill that has gone into the project. ”

Dr Jonathan Hall

Free Site Assessment

To arrange your FREE site assessment, please feel free to call us on: **01322 292222** or email: info@eclipse-dental.com

Guaranteed solutions...

Call us on: **01322 292222** email: info@eclipse-dental.com

www.eclipse-dental.com

Fit-out
Equipment
Service



Ancar dental treatment centres

Comfort and technology, all in one

Irresistible comfort, as well as the easy accessibility of an S7, come together with the great benefits of Ancar's electronic equipment. Enjoy the advantages of the most advanced medical technology, which we call Touch Expert, Ancar's effort to offer the greatest professional services, maximum precision and informative details that come close to perfection.

SD-730 FEATURES

- Whip arm system.
- Vertical elevation chair.
- Joystick seat movements
- Joystick with 4 memories.
- Footrest reclinable to 90° with automatic extension.
- Selection of elevating speeds.
- Electronic equipment, which is suited to the doctor's work and helps to reduce the treatment times thanks to Touch Expert technology.
- 70° movement from the water group to facilitate 4-handed work.
- Junction box built into the side.
- Multifunctional pedal.
- Cannula support arm with adjustable height.
- Touch Expert Technology.



Model shown: SD-730 dental treatment unit

Ancar the most robust chair in the market

- Simplicity and quality of engineering
- Made in Aluminium



Call
01322 292222
now to discuss the
Ancar range of
treatment units

Guaranteed solutions...

Call us on: **01322 292222** email: info@eclipse-dental.com

www.eclipse-dental.com



Space saving air compressor and suction pump

MGF ASPIR-Comp

Featuring:

- Space-saving slim profile
- Quiet operation
- Clean and compact
- Suitable for 1 - 5 surgeries
- Built-in amalgam separation
- Membrane air-dryer
- Oil-free

Contact us today for a brochure.

Call: **01322 292222**

Email: info@eclipse-dental.com



Brighten up your surgery

Zenith LED Light Fittings

Brighten up your surgery with these energy saving LED lights and never have to change bulbs!

- High transmission optics
- Curved reflector
- Superior luminaire in its class
- 4000K colour temperature
- Dimmable option
- Emergency option

Promotional price

£239.00 each +VAT

Installation service available.

Call: **01322 292222**

Email: info@eclipse-dental.com



Guaranteed solutions...

Call us on: **01322 292222** email: info@eclipse-dental.com

www.eclipse-dental.com

Fit-out
Equipment
Service



CrossTex DentaPure® DP365B

Independent water bottle cartridge



- Installs in minutes
- Cost-effective
- Can use tap water
- Works for 365 days

Special offer
Limited stock only

Price:
£149.00 +VAT

Order now

Call: **01322 292222**

Email: info@eclipse-dental.com

One bur does it all



TriHawk Talon metal cutting bur

- Remove 6 crowns in less than 1 minute
- Always cuts regardless of angle
- High cutting speed
- Suitable for removing old crowns, bridges and amalgam
- Cuts through porcelain, metal, zirconium, amalgam and enamel

Pack 25
£2.75
+VAT each

Pack 100
£3.75
+VAT each

Ask us about a sample or demonstration

Call: **01322 292222** Email: info@eclipse-dental.com

EclipseCare

Equipment annual servicing

Our experienced engineers carry out a full, in-depth service fitting manufacturer's recommended service kits to minimise the possibility of break down.

We can even connect a back-up compressor or suction pump during the service to keep your surgery working!



- Simple pricing per item
- PAT testing dental equipment
- Service certificate for CQC inspection
- Planned Preventative Maintenance
- Free, unlimited technical support is available

Call
01322 292222
now to get your
service quote
and schedule

CELEBRATING

12
YEARS

In 2018 Eclipse Dental will be celebrating twelve years, providing dental professionals with guaranteed solutions for...

- Fit-out
- Servicing
- Equipment supplies

If you have any questions about the services we offer, then please feel free to contact us.

Fit-out
Equipment
Service



Eclipse Dental Engineering
Westview Court, 10-14 Kent Road,
Dartford, Kent.DA1 2AJ

Guaranteed solutions...

Call us on: **01322 292222** email: info@eclipse-dental.com

www.eclipse-dental.com